


Fit Your Needs, **GM85 Fit**



FOCUS ON THE CORE OF "MOBILE" RADIOGRAPHY

In recent extraordinary situations, radiography is facing a new paradigm in its ways of working. As the procedure volume of mobile x-ray increases to deal with unprecedented patient volume, Samsung is focusing more on efficient mobile imaging based on its accumulated technology. Experience higher productivity and performance with GM85 Fit designed to fit your workflow, diagnosis and care.



Rapid Charging¹

Within
3~4
hours

All Day Operation²

500
exposures
50 km
(31 miles)
10 hours

Max. Exposure³

2,300
exposures

Sleep Mode Stand-by⁴

31
hours

1) 200 ~ 240 VAC < 3 hours, 100 ~ 127 VAC < 4 hours

2) Test condition : Chest AP / 80 kVp / 160 mA / 5 msec / 60 sec intervals, with driving at 5.6 km/h

3) Test condition : Chest AP / 80 kVp / 160 mA / 5 msec / 10 sec intervals, without driving

4) Stand-by 18 hours for LCD on, 31 hours for sleep mode, and 75 hours for power off

Fit Your Workflow

WHEN YOU GO TO START AN EXAM,



From the moment you start working, GM85 Fit increases your work efficiency with the fast speed of booting up the system. **Quick Driving** allows technologists to move GM85 fit right away without waiting for the system to be fully ready. Users can easily log in to the system through the **embedded RFID badge tag***. It is a built-in type in the equipment, so it reduces external interference.

Additional Patient Information* is useful when checking the list of patients requiring x-ray imaging. It intuitively shows infection and fall risk information and orders comments that support correct action for each patient. This helps keep patients and staff safe by responding to possible risks and taking preemptive measures. The **lightweight** GM85 Fit allows easy maneuvering both in and out of elevators without worrying about the weight limit.

Embedded RFID badge tag*



Additional Patient Information*



Fit Your Workflow

WHEN YOU REACH A PATIENT,




When checking the patient's identity after arriving at the room, a **wireless barcode scanner*** allows you to freely approach the patient and scan information without being restricted by obstacles that may be near you and the patient. With a **wireless hand switch** that supports up to 33 ft (10 m) range, users can take images from far away to reduce radiation exposure.


Since the role of mobile DR is usually to use follow-up exams, it is important to perform x-ray exams under the same conditions. **Prior Exam Review*** displays previous images and exposure parameters of the patient being examined. With a quick comparison, users can improve image consistency and reduce retakes.

SimGrid™* application streamlines the workflow by improving image quality without the use of a conventional grid. This allows the omission of grid installation and removal steps from the conventional workflow leading to 28% reduction in total exam time.

Prior Exam Review*



SimGrid™*



WHEN YOU GET X-RAY IMAGES,



In portable exams, radiologic technologists rotate images after the acquisition as the image orientation depends on the exam environment. **Image Auto rotation***¹ detects the rotated angle of the image and automatically rotates in the correct direction based on AI algorithm. (97% average accuracy) The **auto-run option** of advanced applications creates the companion image without additional settings or x-ray exposure.

CXR Assist* with Lunit INSIGHT CXR Triage prioritizes patients with target indications (Pleural Effusion & Pneumothorax), improving work efficiency and enabling timely care in the ER.

In the ER/OR or Trauma environments, multiple staff may often need to quickly acquire and check images in real-time. **Mirror View*** provides secured screen sharing using Wi-Fi CERTIFIED Miracast. **QuickLink**, you can check whether the acquired images uploaded well to the server by accessing RIS/HIS directly without a separate PC.

When preparing for the next exam, **EasyClean**, which locks the screen for a while, helps you clean equipment quickly and control infection effectively.

CXR Assist*

with Lunit INSIGHT CXR Triage

AI Results		ID	Age	Time
Pef	Ptx			
●		18291CCEE	18	10:50 AM
	●	23075CCEE	63	12:15 PM
		29110CCEE	26	10:12 AM
		31291CCEE	33	10:21 AM
		11007CCEE	43	10:36 AM

This PACS workload is intended to help healthcare professionals understand the functions of this product. The workload may differ from the actual product.


Detector rest

for easy bagging



1) [Chest/Abdomen/Pelvis/Infantogram] 0°, 90°, 180°, 270° Rotation [Hand/Foot] 1° Unit Rotation
2) When using Mirror View, the image transmitted via Miracast is not available for the purpose of diagnosis.

Fit Your Diagnosis




Up to
50%

Dose reduction
Same image quality

Samsung's S-Vue™ image processing engine provides excellent, high resolution images for diagnostic accuracy and confidence. The S-Vue™'s adaptive filtering and advanced denoising technology offer the natural presentation of bone and soft tissues. Also, dose level can be reduced up to 50% for adult chest and 47.5% for adult abdomen without compromise in image quality.

Note : The claim concerning Samsung DR is based on limited phantom and clinical study results. Only routine PA chest radiography and abdominal radiography for average adults and pediatric abdominal, chest, skull radiography were studied, excluding pediatric patients under 1 month old. (FDA cleared - K172229, K180543, K182183) In practice, the values of dose reduction may vary accordingly. These clinical images calculates the dose reduction rate from its own standard dose at the clinical site, unlike our FDA claim which compares dose between new IPE and old IPE. The clinical site is responsible for determining whether the particular radiographic imaging needs are not impacted by such x-ray dose reduction.

Case 1. Chest AP (X-Large adult)



The image was taken with GM85.

Case 2. Abdomen



The image was taken with GC85A.

Case 3. L-spine



The image was taken with GC85A.



Standard Detector

IP57 



S4343-AWM



S4335-AWM



881 lbs (400 kg)
Surface load



441 lbs (200 kg)
Point load¹



441 lbs (200 kg)
Bend pressure

1) Based on 4cm radius on the center

High quality Small pixel pitch

To offer optimized image quality, Samsung introduce AccE Standard Detector that broaden customers' choice of Csi detectors. This latest lineup will meet various clinical needs for sizes.

Enhanced load allowance along with dust and water resistance allows the detector to be actively implemented in versatile environments such as ER and OR.

The new feature, detector built-in AEC, S-AEC enables consistent image quality in mobile imaging by precisely controlling radiation dose according to patient size. It is a truly wireless solution that requires no additional receptors or cables, and provides stable operation without concerns about detector direction, SID, or network latency.

Case 4. Pelvis




The image was taken with GC85A.

Consistent Image Quality

The upgraded S-Vue™ displays areas with varying bone density clearly regardless of the patient's size or body position through advanced region-specialized processing. With this engine, halo artifacts surrounding metallic objects are mitigated and visualized properly.

Fit Your Diagnosis



Auto Lung Nodule Detection^{*1}

Auto Lung Nodule Detection is computer aided detection software to identify and mark regions in relation to suspected pulmonary nodules from 10 to 30 mm in size. It is designed to aid the physician to review the PA chest radiographs of adults as a second reader and be used as part of S-Station.

Case. Chest PA

Images were taken with GC85A.




Without Auto Lung Nodule Detection



With Auto Lung Nodule Detection

¹) ALND cannot be used on the patients who have lung lesions other than abnormal nodules and was not tested on images having more than three nodules.



SimGrid^{TM*}

With just a click, SimGridTM allows you to provide better patient care with higher satisfaction and reduced retake rates without the use of a portable grid. It improves image contrast by reducing scatter radiation effects and creates better image quality. The 3-step intensity control (Low/Medium/High) enables customized image processing.

Case. Chest AP

Images were taken with GM85.



Without Grid



With SimGridTM



S-Enhance*

To support your diagnosis, S-Enhance improves the clarity of foreign bodies (e.g. tube, line and/or needle) in images of chest, abdomen, and L-spine. With a single on-screen click, the companion image is created without additional settings or x-ray exposure, streamlining the workflow.

Case. Chest AP


Images were taken with GM85.



Without S-Enhance



With S-Enhance



Bone Suppression*


Without additional setting or exposure, Bone Suppression Imaging improves the clarity of soft tissues by suppressing the appearance of bones in chest images, which improves your ability to detect nodules. You can easily create the companion image with just a click on the screen.

Case. Chest PA

Images were taken with GC85A.



Without Bone Suppression



With Bone Suppression

Fit Your Care




Up to
45%

Dose reduction For very young patients

Our youngest patients are more radiosensitive than adults. To alleviate these concerns, the new S-Vue™ engine helps achieve the optimal dose level for children during pediatric x-ray scans. The dose level can be reduced up to 45% for pediatric abdomen, 15.5% for pediatric chest, and 27% for pediatric skull exams. This is especially significant as abdomen protocols may include genital regions.

Case. Pediatric Chest AP


Images were taken with GM85.



Conventional
16.1 uGy

(54 kVp / 1.42 mAs / 0.06 dGy*cm² / 0.1 mmCu Filter)

15%
DOSE
REDUCTION




Low Dose
13.7 uGy

(54 kVp / 1.21 mAs / 0.05 dGy*cm² / 0.1 mmCu Filter)

6stages

Pediatric Exposure Management

Optimized **6-stage weight dependent imaging** enables pediatric patients to avoid unnecessary x-ray exposure using precise dose management, resulting in superior image quality. Users who need detailed control can apply 3 types of **additional filters**. (0.1/0.2/0.3 mm Cu)



Less than
43dB

when the system is on stand-by

GM85 Fit has low noise operations that do not disturb other patients. With **night mode** on, the screen color, temperature, and sound volume are accordingly adjusted. It will be a necessary function for taking x-ray exams in a quiet environment such as NICU. Its **quick exposure** feature is useful for pediatric patients who are not easy to maintain stationary posture.



American Academy of Pediatrics (AAP) recommendation in NICU ≤ 45 dB



Efforts to reduce fatigue of radiographers

The **lightweight** GM85 Fit provides an effortless driving experience for users. Not only the equipment but also the detector that must be held directly with the arm should be light. The user-centric design of AccE Standard Detector supports patient positioning and alleviates daily burdens.



747 lb
(339 kg)

Light weight



CENTER ENGRAVING
to help position the patient



REAR GRIP
to support transportation



SIDE CHAMFER
to ease your lifting




Cybersecurity & Value Care Service

Samsung provides a solution to support our customers by offering the tools to protect against cyber threats that may compromise invaluable patient data and ultimately degrade the quality of care. In addition, **remote software update** function allows the system to keep up to date. Experience Samsung's **value care service** that improves the system's uptime and total cost of ownership with 26% faster repair service.¹


Samsung's Defense-in-depth

Designed holistic protection for your critical assets with multi-layered security controls



26%

Faster repair service¹



¹) Based on the service time evaluation for major parts repair (compared to the conventional GM85)

About Samsung Electronics Co., Ltd.

Samsung Electronics Co., Ltd. inspires the world and shapes the future with transformative ideas and technologies.

The company is redefining the worlds of TVs, smartphones, wearable devices, tablets, cameras, digital appliances, medical equipment, network systems, and semiconductor and LED solutions. For the latest news, please visit the Samsung Newsroom at news.samsung.com.

GM85 Fit is a sub-configuration of GM85. (Fit-type)



Scan code or visit
www.samsunghealthcare.com
to learn more

Copyright © 2025 Samsung Electronics Co. Ltd. All rights reserved. Samsung is a registered trademark of Samsung Electronics Co. Ltd. Specifications and designs are subject to change without notice. Non-metric weights and measurements are approximate. All data were deemed correct at time of creation. Samsung is not liable for errors or omissions. All brand, product, service names and logos are trademarks and/or registered trademarks of their respective owners and are hereby recognized and acknowledged. Samsung Electronics Co., Ltd.